



DANASWAB VIRAL™ Sample Collection Kit COVID-19


DANASWAB VIRAL Sample Collection Kit is designed for collection, ambient storage and transport viral DNA/RNA from samples collected using a swab.

Accommodates a wide variety of samples:
Wherever a swab may be deemed appropriate per application.

- *Nasal
- *Throat
- *Saliva
- *Fecal
- *Surfaces
- * Environmental samples

Technology overview:

DANAVIRAL Stabilization Solution effectively inactivates viruses and prevents nucleic acid degradation in samples collected with the provided swab, resulting in non-infectious samples to be handled and shipped safely.



Features

- Inactivate microorganisms and viruses for safe and easy transport and handling.
- Each collection tube is pre-filled with 1 ml.
- Preserve total DNA, including viral DNA, at room temperature (4°-25°C) 1 year.
- Preserve total RNA, including viral RNA, at room temperature (4°-25°C) 1 month.

Ref.0603.47	100 Preps	Stock
Collection tubes with DANAVIRAL Stabilization Solution	1ml x 100 units	Room temperature
Swabs	100 units	Room temperature